

## QT633

### Bed Exit Mat



## Overview

Designed to trigger a call by sensing the pressure change when a patient exits a bed.

Connects to Quantec via the remote jack socket on a QT637 Enuresis/General Purpose Interface Socket or QT636 Interface Unit.

To reduce the risk of nuisance calls, the QT637 includes circuitry that prevents a momentary change in pressure such as a patient quickly changing position from triggering a call.

Typically positioned across and over a mattress, directly underneath the bedsheet.

Measures 80cm x 10cm and is supplied with a 2m lead and normally closed right angled jack plug.

## Technical Specifications

<b>Compatibility</b>	Must be used in conjunction with a QT637 Quantec Enuresis Sensor/General Purpose Interface Socket or QT636 Interface Socket.
<b>Connections</b>	6.35mm (quarter inch) jack plug.
<b>Product dimensions (mm)</b>	766mm W x 102mm H. Includes a 2m lead.
<b>Construction &amp; finish</b>	Plastic.



<b>IP Rating</b>	TBC.
<b>Weight</b>	240g
<b>Operating conditions/temperature</b>	-5°C to +40°C. Max relative humidity: 95%.



T 01942 322744 E sales@c-tec.co.uk  
F 01942 829867 W www.c-tec.com



C-TEC, Challenge Way, Martland Park, Wigan, WN5 0LD United Kingdom



T 01942 322744 E sales@c-tec.co.uk  
F 01942 829867 W www.c-tec.com



C-TEC, Challenge Way, Martland Park, Wigan, WN5 0LD United Kingdom